

Field Semen Evaluation Protocol

Step-by-step checklist for bull semen quality at the point of use



STEP 1 PRE-CHECK — EQUIPMENT & TEMPERATURE

- | | |
|--|--|
| ☐ Thaw unit filled, thermometer in place | ☐ Water bath at 95°F (35—37°C) — confirmed with thermometer |
| ☐ AI gun warmed (rub with paper towel) | ☐ Sheath, scissors, towel, gloves laid out and dry |
| ☐ Cow in chute, calm, ready | ☐ Tank straw located below the frost line |
| ☐ Tweezers ready (never use fingers to pull straws) | |

Note — have everything ready before you thaw. The clock starts the moment the straw leaves the tank.

STEP 2 VISUAL ASSESSMENT

- | | |
|---|--|
| ☐ Straw pulled below frost line in under 5 sec | ☐ Shaken to clear nitrogen from the cotton plug |
| ☐ One straw at a time, 95°F for 30—45 sec | ☐ Dried completely after bath (water is lethal) |
| ☐ Checked for cracks or leaks | ☐ Protected from cold shock — gun at body temp |
| ☐ Inseminated within 6—15 min of thawing | |

Note — never thaw more straws than you can use in 15 minutes. Quality drops by the minute.

STEP 3 MAKSA EVALUATION — OBJECTIVE QUALITY CHECK

- | | |
|---|---|
| ☐ Drop of thawed semen placed on slide | ☐ Total motility recorded (%) |
| ☐ Progressive motility recorded (%) | ☐ Reading meets threshold (min. 30% progressive) |
| ☐ Timestamp + GPS captured automatically | ☐ Record saved — bull ID + batch |
| ☐ If below threshold: do NOT breed — flag the batch and the bull | |

Note — check the dose that reached YOUR farm, not the one tested at the stud weeks ago.

FIELD RECORD

Date	_____	Bull ID	_____
Batch / Straw	_____	Technician	_____
Total motility %	_____	Progressive %	_____
Decision	☐ Breed ☐ Reject		

Want objective motility in the field in 90 seconds? MAKSA turns a smartphone into a semen analyzer — validated at $r > 0.90$ against laboratory CASA systems, used in 10+ countries.

See it work: avareit.com • info@avareit.com • WhatsApp +971 50 641 2775

Based on standard field AI handling practice and Society for Theriogenology semen evaluation guidelines.